

cDNA Array CGH

1. **Purify genomic DNA** (Trizol and phenol/chloroform method)

2. **Sonication**

Using probe sonicator: 10 sec x 6

(20 ug DNA in 300 µl H₂O)

Check on gel--- smear 200-800bp

Purification with Qiagen PCR spin column protocol.

- a. Add 900 µl PB buffer to reaction and transfer to QIAquick column.
- b. Spin at maximum speed for 1 minute.
- c. Wash with 700 µl **PE** buffer
- d. Spin at maximum speed for 1 minute.
- e. Repeat wash steps c and d.
- f. Empty collection tube and spin at maximum speed for 1 minute.
- g. Transfer column to a fresh tube, Elute with 30 µl **H₂O** with 1 min incubation followed by spin at maximum speed for 1 minute.
- h. Repeat elution step g

3. Genomic DNA labeling

For one labeling use 3 μ g of the digested control DNA and test DNA.

Add

DNA 3 μ g	19 μ l
Random Primer (hexamers) Gibco (3 ug/ul)	2 μ l

Incubate for 10 min at 95° C, cool on ice

aadNTPs 50x*	0.5 μ l
Klenow Reaction Buffer 10x	2.5 μ l
Klenow Fragment NEB (50U/μl) **	1 μ l
total	25 μ l

*

dATP 25mM
dCTP 25mM
dGTP 25mM
dTTP 5mM
aa-dUTP 20mM

Note that the ratio between dTTP and aa-dUTP is 1:4

Incubate 37 ° C overnight, then stop reaction by adding 5 μ l 0.5 M EDTA pH 8.0.

Probe purification with Qiagen PCR spin column protocol.

- i. Add 300 μ l PB buffer to reaction and transfer to QIAquick column.
- j. Spin at maximum speed for 1 minute.
- k. Wash with 700 μ l **phosphate wash buffer** (5 mM KPO₄, pH 8.0 + 80% EtOH – use this home-made buffer to avoid Tris).
- l. Spin at maximum speed for 1 minute.
- m. Repeat wash steps c and d.
- n. Empty collection tube and spin at maximum speed for 1 minute.
- o. Transfer column to a fresh tube, Elute with 30 μ l **phosphate elution buffer** (4 mM KPO₄, pH 8.5) with 1 minute incubation followed by spin at maximum speed for 1 minute.
- p. Repeat elution step g.
- q. Speedvac dry (~40 minutes)

4. Cy Dye coupling

- a. Resuspend the cDNA in 4.5 μ l 0.1 M carbonate buffer, pH 9.0 (freshly made, <1 month).
- b. Add 4.5 μ l NHS-Cy dyes (resuspended in DMSO. AVOID ANY MOISTURE!!)
- c. Incubate the reactions in dark for 1 hour in room temperature.

Probe purification with Qiagen PCR spin column protocol.

- a. Add 35 μ l **100 mM** NaAc pH5.2 to the reaction.
- b. Add 300 μ l PB buffer to reaction and transfer to QIAquick column.
- c. Spin at maximum speed for 1 minute.
- d. Wash with 700 μ l **phosphate wash buffer** (5 mM KPO₄, pH 8.0 + 80% EtOH – use this home-made buffer to avoid Tris).
- e. Spin at maximum speed for 1 minute.
- f. Repeat wash steps c and d.
- g. Empty collection tube and spin at maximum speed for 1 minute.
- h. Transfer column to a fresh tube, Elute with 30 μ l H₂O with 1 minute incubation followed by spin at maximum speed for 1 minute.
- i. Repeat elution step g.

At this step, 1 μ l of sample can be taken for dye incorporation analysis. Diluted sample at 1:100. Measure absorbance at 260 nm, 550 nm (for Cy3) and 650 nm (for Cy5).

For Cy3, incorporation: nucleotide/dye=17.1*OD₂₆₀/OD₅₅₀.

For Cy5, incorporation: nucleotide/dye=28.5*OD₂₆₀/OD₆₅₀.

Prehybridization:

- a. Prepare prehybridization buffer (**5xSSC, 0.1% SDS and 1% bovine serum albumin, BSA Sigma Cat# A-9418**)
- b. Place slides to be analysed into a staining dish, fill with prehybridization buffer, and incubate at 42 °C for 45 minutes
- c. Wash the slides in a Wheaton slide rack at room temperature with DIH₂O with agitation.
- d. Dip the slides in room temperature isopropanol and dry by centrifugation at 2000 rpm for 2 min. Inspect slides to make sure there is no protein marks on slides. If there is, repeat step 3 and 4.

Slides should be used immediately following prehybridization, at least within one hour.

Hybridization:

1. Combine the Cy3 and Cy5 labeled probes (about 120 μ l)

add 50 μ g of Cot-1 DNA (10mg/ml) 5 μ l

100 μ g of yeast tRNA (4 mg/ml) 25 μ l

speedvac dry it.

2. Redissolve the combined probes in 36 μ l 1X hybridization solution

Dextran hybridization solution

50% Formamide 500 μ l of Formamide / 1ml

4xSSC 200 μ l of 20xSSC / 1ml

2% SDS 100 μ l of 20% SDS / 1ml

10% Dextran Sulfate 200 μ l of 50% Dex Sul /1ml

Cy5/Cy3 labeled probes	37 μ l
------------------------	------------

3. Denature at 75 °C for 10 minutes

4. Incubate at 37 °C for 1 h.

5. Apply 3x7 μ l H₂O into each modified hybridization chamber to prevent drying.

6. Apply probes to the prehybridized microarray slide, hybridize in a humidified chamber at 37 °C for 16-24 hours.

Wash:

Place slide into a staining dish with Wash #1 until coverslip falls off.

Wash #1:	at R/T in 0.2% SDS + 1 X SSC until coverslip falls off
Wash #2:	4 minutes at 42 °C in 0.2% SDS + 1 X SSC
Wash #3:	4 minutes at R/T in 0.2% SDS + 0.1 X SSC
Wash #4:	4 minutes at R/T in 0.06 X SSC

Spin immediately (1000 rpm for 2 minutes at R/T)

Appendix:

Solutions:

100 mM aa-dUTP (Sigma, cat # A0410)

Disolve 10 mg in 191 ul 0.1 M KPO₄ buffer, pH 7.5. Quantitate this stock solution by diluting an aliquot 1:5000 in the same buffer and measuring A₂₈₉. Stock concentration in mM = OD_{289nm} x 704.

Make 10 ml 0.1 M KPO₄ buffer, pH7.5

1. Prepare 2 solutions: 1 M KH₂PO₄ and 1M K₂HPO₄.

2. Mix:

0.802 ml 1M K₂HPO₄

0.198 ml 1 M KH₂PO₄

3. Bring to 10 ml with H₂O. The pH of this KPO₄ buffer should be at 7.5.

50X aa-dNTP mix

Nucleotide/dye ratio with the ratio of aa-dUTP to dTTP. We use a 2:3 ratio. Prep are a 50X stock:

Start Concentration	Volume	Final Concentration
100 mM dATP	25 µl	25 mM dATP
100 mM dCTP	25 µl	25 mM dCTP
100 mM dGTP	25 µl	25 mM dGTP
100 mM dTTP	5 µl	5 mM dTTP
100 mM dUTP	20 µl	20 mM aa-dUTP

Phosphate wash buffer

1. Prepare 2 solutions: 1 M KH₂PO₄ and 1M K₂HPO₄.

2. Make 1 M KPO₄ (the pH should be 8.5-8.7) mix:

9.5 ml 1M K₂HPO₄

0.5 ml 1 M KH₂PO₄

3. Make 100 ml phosphate washing buffer, Mix:

0.5 ml 1 M KPO₄, pH 8.5

15.25 ml H₂O

84.25 ml 95% EtOH

(This solution will be slightlt cloudy)

Phosphate elution buffer

Dilute 1 M KPO₄ pH 8.5 to 4 mM for Phosphate elution buffer:

1 M KPO ₄ pH 8.5	0.4 ml
H ₂ O	99.6ml
<hr/>	
Total	100ml

Carbonate buffer 0.1 M Na₂CO₃, pH 9.0:

Na ₂ CO ₃	0.27g
H ₂ O	20ml
Adjust pH to 9.0 with 6 N HCl	
Total with H ₂ O	25 ml

NHS-Cy dyes

NHS-Cy3: AmershamPharmacia Cat# PA23001

NHS-Cy5: AmershamPharmacia Cat#PA25001

1. Resuspend one tube in 73 ul DMSO.
2. Use immediately, or aliquot 4.5 ul into 0.5 ml tubes and store at -80°C. Avoid moistures!!!