

Immunocytochemistry - Centrosomal Staining

Section of Cancer Genomics, Genetics Branch, NCI
National Institutes of Health

Reagents

1X Phosphate Buffered Saline (PBS)

GIBCO – Life Technologies, Cat. 10010-023

Bovine Serum Albumin Fraction V, protease free (BSA)

Roche Diagnostics, Cat. 03 117 332 001

4',6-Diamidino-2-phenylindole dihydrochloride (DAPI)

Sigma-Aldrich, Cat. D9542

Ethylene glycol-bis(2-aminoethylether)-*N,N,N',N'*-tetraacetic acid (EGTA)

Sigma-Aldrich, Cat. E4378

1M HEPES Buffer

Sigma-Aldrich, Cat. H3537

2M Magnesium Chloride (MgCl₂)

Quality Biological, Inc., Cat. 340-034-060

Methyl Alcohol (MeOH), Anhydrous (Absolute)

MACRON Fine Chemicals, Cat. 3016-02

Monoclonal Anti-Gamma Tubulin (Mouse)

Sigma-Aldrich, Cat. T 6557

Normal Goat Serum (NGS), Sterile

Sigma-Aldrich, Cat. G6767

PIPES BUFFER (PIPES) No. P6757

Sigma-Aldrich, Cat. P8658

Sheep Anti-Mouse IgG (whole molecule) FITC

Sigma-Aldrich, Cat. F3008

Triton X-100

Sigma-Aldrich, Cat. T8532

dH₂O, sterile

Preparation

PHEM Buffer

PIPES	1.81	g	f.c. [60 mM]
HEPES	0.06	g	f.c. [25 mM]
0.5M EGTA	20	ml	f.c. [10 mM]
2M MgCl ₂	1	ml	f.c. [2 mM]
dH ₂ O, sterile	978	ml	
Total	1000	ml	

*pH to 6.9 with approximately 900 μ l 10M NaOH

0.5 M EGTA

EGTA 3.804 g

Water, sterile 20 ml

*pH to 7.5 to get EGTA into solution

Procedure

1. Wash cells 2x with 1X PBS.
2. Carefully permeabilize cells with 0.5% Triton X-100/PHEM buffer 5 min at room temperature (RT).
3. Carefully wash cells 2x with PHEM.
4. Carefully fix cells with ice-cold (-20°C) MeOH 10 min at RT.
5. Remove sides, if using chamber slides.
6. Wash 4x with 1X PBS shaking for 5 min at RT (for cells attached to chamber slides, wash in coplin jar).
7. Incubate with 200 μ l Blocking solution (5% NGS, 1% BSA in 1X PBS) 30 min at 37°C in moist chamber.
8. Rinse briefly in 1X PBS, add 150 μ l (mouse anti-gamma tubulin diluted 1:1000 in 1% NGS, 1% BSA, 1X PBS) and incubate 45 min at 37°C in moist chamber.
9. Wash 3x 1X PBS 5 min at RT.
10. Add 150 μ l (sheep anti-mouse-FITC diluted 1:200 in 1% NGS, 1% BSA, 1X PBS) and incubate 45 min at 37°C in moist chamber.
11. Wash 3x 1X PBS 5 min at RT.
12. Counterstain DNA with DAPI [80 ng/ml in 2X SSC] 1-5 min at RT.
13. Wash 5-10 min 1X PBS shaking.
14. Mount with antifade and cover slip (invert the cells grown onto coverslips onto a clean glass slide).

Notes

1. Keep chamber slides intact when cells are grown on them throughout Steps 1- 4. If cells were grown on sterile coverslips, these should remain in tissue culture dishes/wells for Steps 1- 13.
2. Steps 1 - 3 can be replaced in some instances by one wash with 1X PBS followed by one wash with PHEM. This is important when doing α -tubulin staining because α -tubulin coagulates into artifactual thick fibers upon permeabilization.
3. This protocol is written for adherent cells grown in chamber slides, but can also be used for fixed cells dropped onto slides (do NOT use hypotonic treated cells if proteins of interest are components of the cytoplasm, as they will be destroyed upon dropping cells onto slide). Slides should be allowed to dry completely and then fixed directly in ice-cold MeOH (-20°C) for 10 min (step 4).