

Immunocytochemistry Followed by Electron Microscopy

Section of Cancer Genomics, Genetics Branch, NCI
National Institutes of Health

Reagents

Antifade (1.4-phenylene-diamine)

Bovine Serum Albumin (BSA)

Boehringer Mannheim Biochemicals (BMB), Cat. 100 350

Gridded coverslips

Bellco, Cat. 1916-92525

DAPI

BMB, Cat. 236 276

EGTA

Sigma, Cat. E3889

Glutaraldehyde, 25% EM grade

Polysciences, Inc., Cat. 01909

2° Antibodies (Ab)

Specific for 1° antibodies

Normal Goat Serum (NGS)

Sigma, Cat. G6767

Magnesium chloride (MgCl₂), 2 M

Quality Biological, Inc., Cat. 340-034-060

1X Phosphate Buffered Saline (PBS), pH 7.4

Invitrogen Corp., Cat. 10010-023

PIPES

Sigma, Cat. P9291

1° Antibodies (Ab)

Specific for desired protein

Sodium borohydride (NaBH₄)

Sigma, Cat. S9125

Sodium chloride (NaCl)

Mallinckrodt, Cat. 7581

Sucrose

Sigma, Cat. S7903

Triton X-100

Calbiochem, Cat. 648462

Preparation

Cytoskeleton (CSK) Buffer

PIPES	1.512 g	f.c. [10mM]
Sucrose	51.35 g	f.c. [300mM]
NaCl	2.923 g	f.c. [100mM]
0.5 M EGTA	0.19 g	f.c. [1mM]
2M MgCl ₂	0.75 ml	f.c. [3mM]

*Bring to 500 ml with sterile distilled water

0.1% Sodium Borohydride solution

Prepared fresh 1mg/ml in 1X PBS

Blocking Solution (5% NGS/ 5% BSA/1X PBS)

NGS	500 µl
BSA	0.5 g
1X PBS	10 ml

Store at 4°C

Antibody Solution (1% NGS/1% BSA/1X PBS)

Blocking Solution I	200 µl
1X PBS	800 µl

DAPI (stock solution)

DAPI	2 mg	f.c. [0.2 mg/ml]
dH ₂ O	10 ml	

Aliquot and store at -80°C

DAPI (staining solution)

DAPI stock solution	40 µl	f.c. [80 ng/ml]
2X SSC	100 ml	

Antifade (1,4-phenylene-diamine)

See Antifade preparation procedure in CGH Protocols

Procedure

1. Grow adherent cells on gridded coverslips or cytopsin suspension cells onto poly-L-lysine coated gridded coverslips.
2. Fix cells in CSK Buffer + 0.1% Triton-X100 for 30 sec.
3. Fix cells in CSK Buffer + 0.1% Triton-X100 + 1% glutaraldehyde for 2 min.
4. Fix cells in CSK Buffer + 1% glutaraldehyde for 10 min.
5. Wash 2 x 15 min with fresh 0.1% sodium borohydride solution.
6. Block coverslips with 25 μ l blocking solution in hybridization chamber 30 min at 37°C.
7. Incubate with 1° Ab in 25 μ l antibody solution in hybridization chamber at 37°C for 60 min.
8. Wash 3 x 5 min with 1X PBS at RT.
9. Incubate with 2° Ab in 25 μ l antibody solution in hybridization chamber at 37°C for 60 min.
10. Wash 3 x 5 min with 1X PBS at RT.
11. Stain for 2 min with DAPI.
12. Wash in 1X PBS for 10 min, shaking.
13. Mount coverslip with 10 μ l antifade on microscope slide.
14. Image cells of interest.
15. Send coverslips in 1X PBS for EM analysis.