

## **Immunocytochemistry Followed by FISH (Version 3)**

Section of Cancer Genomics, Genetics Branch, NCI  
National Institutes of Health

**\*We present multiple versions of the Immunocytochemistry / FISH protocols because the conditions under which each antibody and DNA probe work can be different and must be determined empirically.**

### **Reagents**

**Acetic acid, glacial**

**Bovine Serum Albumin (BSA)**

Roche Diagnostics, Cat. 100350

**DAPI**

Sigma-Aldrich, Cat. 18860

**Formamide**

Sigma-Aldrich, Cat. 47670

**Goat anti-rabbit-TRITC (secondary)**

Sigma-Aldrich, Cat. T-5268

**HCl, 1 M**

**Methanol**

**Para-Formaldehyde**

Sigma-Aldrich, Cat. P6148

**Phosphate Buffered Saline, pH 7.4**

Life Technologies, Cat. 10010-023

**Primary antibody**

Specific for desired protein, made in either a mouse or rabbit

**Rabbit anti-mouse-TRITC (secondary)**

Sigma-Aldrich, Cat. T2402

**NaOH, 0.1 M**

**20X SSC**

**Tween 20**

Sigma-Aldrich, Cat. P1379

**Vysis CEP® Probe**

Vysis (Abbott Molecular)

## Preparation

### Methanol

1. Room temperature
2. Pre-chill to -20°C

### Permeabilization Buffer

Triton X-100	50 $\mu$ l	f.c. [0.53]
1X PBS	10 ml	

### Blocking Solution (3% BSA/1X PBS)

BSA	0.3 g
1X PBS	10 ml

Store at 4°C

### Antibody Solution (1% BSA/1X PBS)

Blocking solution	300 $\mu$ l
1X PBS	600 $\mu$ l

### 2% p-formaldehyde

p-formaldehyde	2 g	
1X PBS	100 ml	
0.1 N NaOH	500 $\mu$ l	f.c. [0.5 mM]

Adjust to pH 7.4 with HCl  
Store <1 month at 4°C

### 50% FA/SSC

20X SSC	20 ml
dH <sub>2</sub> O	80 ml
Formamide	100 ml
<hr/> Total	<hr/> 200 ml

Adjust pH to 7-7.5 with 1 M HCl

**Pre-warm to 45°C**

### DAPI (stock solution)

DAPI	2 mg	f.c. [0.2 mg/ml]
dH <sub>2</sub> O	10 ml	

Aliquot and store at -80°C

**DAPI (staining solution)**

DAPI stock solution	40 µl	f.c. [80 ng/ml]
2X SSC	100 ml	

**Antifade** (1,4-phenylene-diamine)

See Antifade preparation procedure in CGH Protocols

**Procedure**

1. Grow adherent cells in chamber slides.
2. Fix cells in methanol pre-chilled to -20°C for 10 min at RT.
3. Wash 3 x 5 min 1X PBS at RT.
4. Permeabilize cells with 0.5% Triton X-100/PBS 5 min at RT.
5. Wash 3 x 5 min 1X PBS at RT.
6. Remove chamber and block slides with 120 µl blocking solution in hybridization chamber 30 min at 37°C.
7. Incubate with 1° Ab (rabbit or mouse) in 120 µl antibody solution in hybridization chamber at 37°C for 45 min.
8. Wash 3 x 5 min with 1X PBS at RT.
9. Incubate with 2° Ab [goat anti-rabbit-TRITC (1:200) or Rabbit anti-mouse-TRITC (1:200), respectively, in 120µl antibody solution] in hybridization chamber at 37°C for 60 min.
10. Wash 3 x 5 min 1X PBS at RT.
11. Fix with methanol:acetic acid (3:1) at RT 10 min.
12. 2% p-formaldehyde at RT for 1 min.
13. 70%, 90%, 100% ethanol series (3 min each).

**Note: Can counterstain with DAPI at this point and mount slides with antifade to have a look at them and determine if Ab detection worked. Wash coverslips 3 x 5 min in 2X SSC before continuing with procedure.**

14. Combine 1 µl Vysis CEP® probe, 1 µl water, and 7 µl Vysis Hyb Buffer.
15. Add probe cocktail to slide, coverslip, and seal with rubber cement.

16. Denature slide 75°C for 5 min on slide warmer.
17. Incubate in hybridization chamber at 37°C overnight.
18. Remove rubber cement.
19. Wash in FA/SSC pre-warmed to 45°C for 21 min, shaking.
20. Stain for 2 min with DAPI.
21. Wash in 2X SSC for 10 min, shaking.
22. Mount with antifade.